

# Aeroocom's 6" Antimicrobial Leak Proof Carrier

## Leak Proof Advantage:

Eliminate infectious control worries with Aeroocom's Anti-microbial Leak Proof Carrier. Keep spills or leakage contained inside the leak proof carrier preventing system contamination.



## IONPURE Additive:

All Aeroocom healthcare leak proof carriers are manufactured with an additive called Ionpure.

Ionpure is an inorganic additive of soluble glass powder containing anti-microbial metal ions.

## Specifications:


OD D = 6" / H = 14.76"  
Payload D = 4.0" / H = 12.4"  
Drive Rings: Velcro  
Weight: 2.97lbs / 1,350kg  
RFID upon request

## IONPURE Additive:


Glass powder composition combined with Anti-microbial metal ions. Durability and effectiveness has no wash deterioration. Heat Transforming Temperature Load 18.5kg is 132.8°C @ .5% Ionpure. Impact Strength is 83.86kg · cm/cm<sup>2</sup>. Transparency maintains its refractive index.

## IONPURE Testing Results after 1 year

Product made with 2% Ionpure additive



Product made without Ionpure additive



## 3 Easy Steps—One Handed Opening:



1. Place fingers in the indented finger port on the left side—push out.
2. Place thumb under clasp on right side—flip up.
3. With fingers in port, push out and carrier top will swivel open.

**aeroocom**

**Leak Proof**

**Anti-Microbial**

**Infectious Control**